

RESULTS INTERPRETATION

The test results should not be read after 30 minutes.

Multiple Respiratory Multipathogen Antigen Test Kit

(Immunochromatographic Assay) ICA2114-7


[POSITIVE]



SARS-COV-2(COVID19)POSITIVE



INFLUENZA A (FLUA) POSITIVE




INFLUENZA B (FLUB) POSITIVE



RESPIRATORY SYNCYTIAL
VIRUS (RSV) POSITIVE




ADENOVIRUS (ADV) POSITIVE



MULTIPLE POSITIVE

[NEGATIVE]

Two dark blue/purple lines appear in the (C) section and no line appears in the detection area.



[INVALID]

If any (C) line do not appear, the test results must be interpreted as invalid.

